# Symposium Potentials for NVs sensing magnetic phases, textures and excitations (SYNV)

jointly organised by the Low Temperature Physics Division (TT), the Magnetism Division (MA), the Metal and Material Physics (MM), and the Semiconductor Division (HL)

Hans Huebl Walther-Meissner-Institut Bayerische Akademie der Wissenschaften Walther-Meißner-Straße 8 D-85748 Garching hans.huebl@tum.de Matthias Althammer Walther-Meissner-Institut Bayerische Akademie der Wissenschaften Walther-Meißner-Straße 8 D-85748 Garching matthias.althammer@wmi.badw.de Andreas Stier Walter Schottky Institut Technische Universität München Am Coulombwall 4 D-85748 Garching andreas.stier@wsi.tum.de

The understanding of ordered magnetic phases is essential from fundamental and application perspectives. While ferromagnets are currently used in non-volatile, high-density data storage devices, topologically trivial and non-trivial magnetic phases are discussed for future low-energy information processing schemes. An overarching scientific quest towards achieving high performance devices, is to obtain a detailed understanding of the stabilization of magnetic phases, the interactions at play, the resulting magnetic textures, and their excitations. Color centers such as nitrogen vacancies (NV) in diamond represent a modern platform to access key magnetization parameters on the nanoscale and therefore represent an important method to provide a deeper insight into the mechanisms at play. This symposium brings together experts from different but overlapping fields of static and dynamic NV sensing, using tip-based single NV scanning methods as well as ensemble-based widefield imaging concepts. The speakers will provide an overview of the recent and exciting developments of this spectroscopy technique with a particular emphasis on material science questions.

## Overview of Invited Talks and Sessions

(Lecture hall Audimax 2)

### Invited Talks

SYNV 1.1	Mon	13:30-14:00	Audimax $2$	Harnessing Nitrogen Vacancy Centers in Diamond for Next-
				Generation Quantum Science and Technology — •CHUNHUI DU
SYNV 1.2	Mon	14:00-14:30	Audimax $2$	Nanoscale imaging of spin textures with single spins in diamond
				— •Patrick Maletinsky
SYNV 1.3	Mon	14:30-15:00	Audimax 2	Spin-based microscopy of 2D magnetic systems - •JÖRG
				WRACHTRUP
SYNV 1.4	Mon	15:15-15:45	Audimax 2	Exploring antiferromagnetic order at the nanoscale with a single
				spin microscope — •Vincent Jacques
SYNV 1.5	Mon	15:45 - 16:15	Audimax 2	Nanoscale magnetic resonance spectroscopy with NV-diamond
				quantum sensors — • Dominik Bucher

#### Sessions

SYNV 1.1–1.5	Mon	13:30-16:15	Audimax 2	Potentials for NVs sensing magnetic phases, textures an	d
				excitations	

Location: Audimax 2

#### SYNV 1: Potentials for NVs sensing magnetic phases, textures and excitations

Time: Monday 13:30-16:15

# Invited TalkSYNV 1.1Mon 13:30Audimax 2Harnessing Nitrogen Vacancy Centers in Diamond for Next-<br/>Generation Quantum Science and Technology — •CHUNHUI DU<br/>— University of California, San Diego, California, United State

Advanced quantum systems are integral to both scientific research and modern technology enabling a wide range of emerging applications. Nitrogen vacancy (NV) centers, optically active atomic defects in diamond, are directly relevant in this context due to their single-spin sensitivity and remarkable functionality over a broad temperature range. In this talk, I will present our recent efforts on developing NV-based quantum sensing and imaging techniques and their potential to address the challenges in both condensed matter physics and quantum science and technologies. Specifically, we have achieved: 1) electrical control of coherent spin rotations of NV spin qubits in NV-magnon based hybrid systems, 2) nanoscale imaging of magnetic flux and magnetization of a topological superconductor by NV wide-field microscopy, and 3) local non-invasive measurements of thermal environment of Mott insulators by NV centers. Our results illustrate the unique capability enabled by NV centers in exploring the exotic spin, heat, and charge transport in emergent quantum materials. The demonstrated coupling between NV centers and magnons in hybrid quantum devices further points to the possibility to establish macroscale entanglement between distant spin qubits and paves the way for developing transformative NV-based quantum computer.

Invited Talk SYNV 1.2 Mon 14:00 Audimax 2 Nanoscale imaging of spin textures with single spins in diamond — •PATRICK MALETINSKY — Department of Physics, Basel University, Switzerland

Quantum two-level systems offer attractive opportunities for sensing and imaging at the nanoscale. Since its inception, this idea [1] has advanced from proof of concept [2] to a mature quantum technology [3], which finds applications in condensed matter physics, materials science and engineering. In this talk, I will discuss our approach to realizing such quantum sensors [4] and highlight some particularly rewarding applications in the imaging of nanoscale spin textures.

Specifically, I will discuss how we employ single electronic spins in diamond for nanoscale probing of antiferromagnetic systems [5,6] and high-resolution imaging of atomically thin van der Waals magnets [7]. For both, the combination of sensitivity, spatial resolution and quantitative imaging enables unprecedented insights into nanoscale domains and domain-walls down to the atomic monolayer limit.

I will conclude with an outlook of future developments of single spin magnetometers for extreme conditions, such as high magnetic fields, mK temperatures or high-frequency sensors to probe the dynamics of nanomagnetic systems.

 B. Chernobrod et al., J. Appl. Phys. 97, 014903 [2] G. Balasubmaranian et al., Nature 455, 644 [3] P. Appel et al., Rev. Sci. Instr. 87, 063703 [4] N. Hedrich et al. Phys. Rev. Appl. 14, 064007 [5] P. Appel et al., Nano Lett. 19, 1682 [6] N. Hedrich et al., Nature Phys. 17, 574 [7] L. Thiel et al., Science 364, 973

Invited TalkSYNV 1.3Mon 14:30Audimax 2Spin-based microscopy of 2D magnetic systems — •JÖRGWRACHTRUP — 3rd Physics Institute, University of Stuttgart,<br/>Stuttgart, Germany

The investigation of magnetic order in 2D materials requires dedicated

probes. While conventional probes of magnetism with nanoscale resolution, like Lorenz microscopy or MFM fail for few layer- or monolayer samples, STM requires dedicated sample preparation. NV-based magnetic probes on the other hand are very well suited to provide quantitative data with a few ten nm spatial resolution and sufficient sensitivity, even for monolayer samples. In the talk I will describe experiments on CrBr3 which show the domain structure of the material [1]. Upon imaging material with different number of layers we gained insight into interlayer coupling and it impact on magnetic order. We also measure magnetic order over a wide range of temperatures and derive information on the physics of the phase transition of CrBr3. In addition, I will show measurements on CrI3 samples of different thickness and relative orientation where we find signatures of Moiré patterns in twisted multilayers [2].

[1] Qi-Chao Sun et al. Magnetic domains and domain wall pinning in atomically this CrBr3 revealed by nanoscale imaging, Nature Comm. 12, 1989 (2021) [2] Qi-Chao Sun et al. Direct visualization of magnetic domains and moiré magnetism in twisted two-dimensional magnets, submitted (2021)

#### $15~\mathrm{min.}$ break

Invited TalkSYNV 1.4Mon 15:15Audimax 2Exploring antiferromagnetic order at the nanoscale with asingle spin microscope•VINCENT JACQUESLaboratoireCharles Coulomb, Université de Montpellier and CNRS, 34095Montpellier, France

Experimental methods allowing for the detection of single spins in the solid-state, which were initially developed for quantum information science, open new avenues for the development of highly sensitive quantum sensors. In that context, the electronic spin of a single nitrogen-vacancy (NV) defect in diamond can be used as an atomicsized magnetometer, providing an unprecedented combination of spatial resolution and magnetic sensitivity under ambient conditions [1]. In this talk, I will illustrate how scanning-NV magnetometry can be used as a powerful tool for exploring condensed-matter physics, focusing on chiral spin textures in antiferromagnetic materials [2,3].

References: [1] L. Rondin et al., Rep. Prog. Phys. 77, 056503 (2014) [2] I. Gross et al., Nature 549, 252 (2017) [3] A. Finco et al., Nat. Comm. 12, 767 (2021)

Invited Talk SYNV 1.5 Mon 15:45 Audimax 2 Nanoscale magnetic resonance spectroscopy with NVdiamond quantum sensors — •DOMINIK BUCHER — Technical University of Munich, Physical Chemistry, Garching, Germany

Recently, optically probed nitrogen-vacancy (NV) point defects in diamond have emerged as a new class of quantum sensors allowing the detection of magnetic fields on unprecedented length scales. This technique allows the measurement of magnetic resonance signals on the nanoscale down to a single electronic or nuclear spin. In the first part of this talk, I will introduce the concept of quantum sensing with NV-centers, in particular, the detection of oscillating magnetic fields, important for magnetic resonance spectroscopy applications. In the second part, I will report on recent progress on probing thin films and 2D materials with this technique. In the concluding part, possible applications in sensing magnetic phases, textures and excitations will be discussed.